

littleSTEPS® Foot Orthotics

CORRECT FIT:



Ends just proximal to the metatarsal heads



Ends just proximal to the metatarsal heads


Note that children's shoe sizes tend to vary quite a bit between shoe vendors. This chart is meant to be a guide to fitting dimensions for finding the appropriate littleSTEPS® foot orthotic size for your child. littleSTEPS® foot orthotic sizing may not directly correlate with the actual shoe size your child is wearing, but should be relatively close.

Sizing Chart

SIZE	MODEL	U.S. KIDS	U.S. MEN	U.S. WOMEN	EURO
9	YOUTH/ADULT*	7.5 to 8.5	7.5 to 8.5	9 to 10	41 to 42
8	YOUTH/ADULT*	6 to 7	6 to 7	7.5 to 8.5	39 to 40
7	YOUTH/ADULT*	4.5 to 5.5	4.5 to 5.5	6 to 7	37 to 38
6	CHILDREN	3 to 4			35 to 36
5	CHILDREN	1.5 to 2.5			33 to 34
4	CHILDREN	13 to 1			31 to 32
3	CHILDREN	11.5 to 12.5			29 to 30
2	CHILDREN	10 to 11			27 to 28
1	TODDLER	8.5 to 9.5			25 to 26
0	TODDLER	7 to 8			23 to 24
00	TODDLER	5 to 6.5			21 to 22

*Youth sizes are roughly equivalent to Men's sizes in length.

littleSTEPS® Gait Plates



CORRECT FIT:





Inside Arch of Left Foot



Outside of Left Foot

Note that children's shoe sizes tend to vary quite a bit between shoe vendors. This chart is meant to be a guide to fitting dimensions for finding the appropriate littleSTEPS® foot orthotic size for your child. littleSTEPS® foot orthotic sizing may not directly correlate with the actual shoe size your child is wearing, but should be relatively close.

Sizing Chart

Toddler

SIZE	MODEL	U.S. KIDS	U.S. MEN	U.S. WOMEN	EURO
7	YOUTH/ADULT*	4.5 to 5.5	4.5 to 5.5	6 to 7	37 to 38
6	CHILDREN	3 to 4			35 to 36
5	CHILDREN	1.5 to 2.5			33 to 34
4	CHILDREN	13 to 1			31 to 32
3	CHILDREN	11.5 to 12.5			29 to 30
2	CHILDREN	10 to 11			27 to 28
1	TODDLER	8.5 to 9.5			25 to 26
0	TODDLER	7 to 8			23 to 24
00	TODDLER	5 to 6.5			21 to 22

*Youth sizes are roughly equivalent to Men's sizes in length.